

MicroFit HealthWizard - Version 5.3.6

May 2010

Fitness Profile Bug Fixes

CTI communication problems with Windows 7 and USB to Serial Adapter

Communication problems with the CTI (white box), blood pressure graphics and aerobic heart rate graphics occurred when using a USB to serial adapter with a PC running Windows 7 OS.

Blood pressure text in Fitness Profile report not showing when BP is 160/100 or above

When systolic blood pressure is 160 mmHg or above or diastolic blood pressure is 100 mmHg or above, the text in the Fitness Profile report for the Needs Work category is not printing.

Default bike type configuration is Monark 818E

When HealthWizard is installed for the first time, the software is configured for the Monark 818e bike. It has been changed to default to the Monark 828e bike.

**Please go to our website to learn the new features from earlier version updates:
http://www.microfit.com/site_fitness/support/updates.html**