

MicroFit HealthWizard - Version 5.3.7

February 17, 2011

Fitness Profile

Bug Fix: Total Fitness Score does not equal 100
When a client achieved the maximum (best) scores for all the fitness tests in a Total Fitness Score calculation, the Total Fitness Score was shown as 99 when it should be 100.

Bug Fix: Edit Aerobic Fitness score requires height
If a Fitness Profile format is created with only the weight and aerobic fitness tests (no height) and a client is tested with this format using the U.S.A.F. Bike protocol, when you attempt to edit the client's scores for this test session the edit program shows a message stating the test requires a height measurement and does not let you edit the scores. With this software upgrade, under this scenario a message box now appears requesting height data entry.

Wellness Profile

Bug Fix: Wellness Profile report has bad exercise statement
When a client answered NO to the question 'Do you warm up before and cool down after aerobic exercise?' the Wellness Profile report stated the client HAD established the important habit of warm-up and cool-down. This was an incorrect statement.

HealthWizard

Bug Fix: Print reports to file: error 75 – file/path access error
When printing any HealthWizard report to a PDF file, there was a 'File/Path Error' and the report was not created. Note: this error occurred only for certain versions of HealthWizard and only if the user's application folder was not set for read/write access.